

Product Datasheet

Human 2B4 / CD244 / SLAMF4 Protein (orb257160)

Description Human 2B4, His Tag is expressed from human 293 cells (HEK293). It contains AA Cys 22 - Arg 221 (Acce...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target 2B4

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 23.1 kDa. The protein migrates as 35-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_057466.1

Purity 95%


Source Human 2B4, His Tag (orb257160) is expressed from human 293 cells (HEK293). It contains AA Cys 22 - Arg 221 (Accession # Q9BZW8-2).

MW 23.1 kDa

Uniprot ID **Q9BZW8**

NCBI **NP_057466.1**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human 2B4
protein