

Product Datasheet

Human 2B4 / CD244 / SLAMF4 Protein (orb257161)

Description Human 2B4, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Cys 22 - Arg 221 (Accession # Q9BZW8-2).

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target 2B4

Preservatives 50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 48.9 kDa. The protein migrates as 60-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_057466.1

Purity 95%


Source Human 2B4, Fc Tag (orb257161) is expressed from human 293 cells (HEK293). It contains AA Cys 22 - Arg 221 (Accession # Q9BZW8-2).

MW 48.9 kDa

Uniprot ID [Q9BZW8](#)

NCBI [NP_057466.1](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human 2B4 protein