

Product Datasheet

Human 4-1BB / TNFRSF9 Protein (orb257164)

Description Human 4-1BB, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 24 - Gln 186 (Acc...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target 4-1BB

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 43.4 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_001552.2

Purity 95%


Source Human 4-1BB, Fc Tag (orb257164) is expressed from human 293 cells (HEK293). It contains AA Leu 24 - Gln 186 (Accession # NP_001552.2).

MW 43.4 kDa

Uniprot ID [Q07011](#)

NCBI [NP_001552.2](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human 4-
1BB protein