

Product Datasheet

Human Activin RIIB / ACVR2B Protein (orb257168)

Description Human Activin RIIB, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 19 - Thr ...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Activin RIIB

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 14.5 kDa. The protein migrates as 20-35 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_001097.2

Purity 90%


Source Human Activin RIIB, His Tag (orb257168) is expressed from human 293 cells (HEK293). It contains AA Ser 19 - Thr 137 (Accession # Q13705-1).

MW 14.5 kDa

Uniprot ID [Q13705](#)

NCBI [NP_001097.2](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
Activin RIIB
...