

Product Datasheet

Human BMPR-1A / ALK-3 Protein (orb257183)

Description

Human BMPR-1A, Fc Tag (AL3-H5259) is expressed

24 ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

BMPR-1A

Preservatives

50 mM Tris, 100 mM Glycine, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 40.3 kDa. The protein migrates as 56 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

CAA80257

Purity

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human BMPR-1A, Fc Tag (orb257183) is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Arg 152 (Accession # P36894-1).

MW

40.3 kDa

Uniprot ID


[P36894](#)

NCBI

[CAA80257](#)

Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human BMPR-1A prote...