


Product Datasheet

Human Angiopoietin-like 3 / ANGPTL3 Protein (orb257186)

Description	Angiopoietin-like protein 3 (ANGPTL3) is also known as Angiopoietin-related protein 3, Angiopoietin-...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Angiopoietin-like 3
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 24.6 kDa. The protein migrates as 28-35 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_055310
Purity	95%
Source	Human Angiopoietin-like 3, His Tag (orb257186) is expressed from human 293 cells (HEK293). It contains AA Ser 17 - Pro 220 (Accession # NP_055310).
MW	24.6 kDa
NCBI	NP_055310
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
Angiopoietin-
...