

Product Datasheet

Human Annexin A5 Protein (orb257188)

Description Human ANXA5, Tag Free (AN5-H5118) is expressed from E. coli cells. It contains AA Met 1 - Asp 320 (A...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Annexin A5

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag Native

Note For research use only

Application notes This protein carries no "tag". The protein has a calculated MW of 35.9 kDa. The protein migrates as 33-35 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence AAH01429

Purity ● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.


Source Human ANXA5, Tag Free (orb257188) is expressed from E. coli cells. It contains AA Met 1 - Asp 320 (Accession # P08758-1).

MW 35.9 kDa

Uniprot ID **P08758**

NCBI **AAH01429**

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human Annexin A5 pr...