


Product Datasheet

Human Annexin A5 Protein (orb257189)

Description	Human Annexin A5, His Tag is expressed from E. coli cells. It contains AA Met 1 - Asp 320 (Accession...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Annexin A5
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 36.8 kDa. The protein migrates as 35 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	AAH01429
Purity	98%
Source	Human Annexin A5, His Tag (orb257189) is expressed from E. coli cells. It contains AA Met 1 - Asp 320 (Accession # P08758-1).
MW	36.8 kDa
Uniprot ID	P08758
NCBI	AAH01429
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
Annexin A5
pr...