

Product Datasheet

Human Axl Protein (orb257196)

Description Human Axl, His Tag is expressed from human 293 cells (HEK293). It contains AA Ala 26 - Pro 449 (Acce...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Axl

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 46.6 kDa. The protein migrates as 60-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAH32229

Purity 95%


Source Human Axl, His Tag (orb257196) is expressed from human 293 cells (HEK293). It contains AA Ala 26 - Pro 449 (Accession # AAH32229).

MW 46.6 kDa

Uniprot ID **P30530**

NCBI **AAH32229**

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human Axl protein