

Product Datasheet

Human Azurocidin / CAP37 / AZU1 Protein (orb257197)

Description

Human Azurocidin, His Tag is expressed from



Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Azurocidin

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-6×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 25.0 kDa. The protein migrates as 34-42 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

AAH69495

Purity

95%

Source

Human Azurocidin, His Tag (orb257197) is expressed from human 293 cells (HEK293). It contains AA Ile 27 - Pro 250 (Accession # AAH69495).

MW

25.0 kDa

Uniprot ID


[P20160](#)

NCBI

[AAH69495](#)

Expiration Date

6 months from date of receipt.



SDS-PAGE
analysis of
Human
Azurocidin
pr...