

Product Datasheet

Human B7-H5 / Gi24 / VISTA Protein (orb257209)

Description Human B7-H5, His Tag is expressed from human 293 cells (HEK293). It contains AA Phe 33 - Ala 194 (Ac...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target B7-H5

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-10×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 20.0 kDa. The protein migrates as 33-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAH20568

Purity 90%


Source Human B7-H5, His Tag (orb257209) is expressed from human 293 cells (HEK293). It contains AA Phe 33 - Ala 194 (Accession # AAH20568).

MW 20.0 kDa

Uniprot ID [Q9H7M9](#)

NCBI [AAH20568](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human B7-H5 protein