


Product Datasheet

Human Bcl-x / BCL2L1 Protein (orb257217)

Description	Human Bcl-x, GST Tag is expressed from E. coli cells. It contains AA Ser 2 - Arg 212 (Accession # AA...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Bcl-x
Preservatives	20 mM Tris, 100 mM NaCl, pH7.5 with 10 mM GSH
Form/Appearance	Powder
Storage	-20°C
Tag	N-GST
Note	For research use only
Application notes	This protein carries a GST tag at the N-terminus. The protein has a calculated MW of 50.6 kDa. The protein migrates as 45-50 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	AAH19307
Purity	92%
Source	Human Bcl-x, GST Tag (orb257217) is expressed from E.coli cells. It contains AA Ser 2 - Arg 212 (Accession # AAH19307).
MW	50.6 kDa
Uniprot ID	Q07817
NCBI	AAH19307
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human Bcl-
x protein