


Product Datasheet

Human CD48 / SLAMF2 Protein (orb257218)

Description	Human CD48, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 27 - Ser 220 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD48
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 23.2 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH16182
Purity	95%
Source	Human CD48, His Tag (orb257218) is expressed from human 293 cells (HEK293). It contains AA Gln 27 - Ser 220 (Accession # AAH16182).
MW	23.2 kDa
Uniprot ID	P09326
NCBI	AAH16182
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
CD48
protein