

Product Datasheet

Human IL-2 R beta / CD122 Protein (orb257259)

Description

Human IL-2 R beta Protein, His Tag is expressed


Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-2 R beta

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-6×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 25.4 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source


Human IL-2 R beta Protein, His Tag (orb257259) is expressed from human 293 cells (HEK293). It contains AA Ala 27 - Asp 239 (Accession # NP_000869.1).

MW


25.4 kDa

Uniprot ID
P14784
NCBI
NP_000869.1
Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human IL-2 R beta p...



Human IL-2 R beta Protein, His Tag (Cat....