

Product Datasheet

Human CD72 Protein (orb257261)

Description Human CD72, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Arg 117 - Asp 359 (Acc...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target CD72

Preservatives 50 mM Tris, 100 mM Glycine, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 54.2 kDa. The protein migrates as 60-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAH30227

Purity 95%


Source Human CD72, Fc Tag (orb257261) is expressed from human 293 cells (HEK293). It contains AA Arg 117 - Asp 359 (Accession # AAH30227).

MW 54.2 kDa

Uniprot ID **P21854**

NCBI **AAH30227**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human CD72
protein