


Product Datasheet

Human CD83 Protein (orb257262)

Description	Human CD83, His Tag is expressed from human 293 cells (HEK293). It contains AA Thr 20 - Ala 143 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD83
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 14.9 kDa. The protein migrates as 17-18 kDa, 20-22 kDa and 25-30 kDa under reducing (R) condition (SDS-PAGE) due to different glycosylation.
Protein Sequence	AAH30830.1
Purity	95%
Source	Human CD83, His Tag (orb257262) is expressed from human 293 cells (HEK293). It contains AA Thr 20 - Ala 143 (Accession # Q01151-1).
MW	14.9 kDa
Uniprot ID	Q01151
NCBI	AAH30830.1
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human CD83
protein