


Product Datasheet

Human CD14 Protein (orb257270)

Description	Human CD14, His Tag is expressed from human 293 cells (HEK293). It contains AA Thr 20 - Met 344 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD14
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 35.9 kDa. The protein migrates as 40-53 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH10507
Purity	98%
Source	Human CD14, His Tag (orb257270) is expressed from human 293 cells (HEK293). It contains AA Thr 20 - Met 344 (Accession # AAH10507).
MW	35.9 kDa
Uniprot ID	P08571
NCBI	AAH10507
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human CD14
protein