

Product Datasheet

Human CD55 / DAF Protein (orb257279)

Description Human CD55, His Tag is expressed from human 293 cells (HEK293). It contains AA Asp 35 - Ser 353 (Acc...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target CD55

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 35.8 kDa. The protein migrates as 50-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_000565.1

Purity 95%


Source Human CD55, His Tag (orb257279) is expressed from human 293 cells (HEK293). It contains AA Asp 35 - Ser 353 (Accession # NP_000565.1).

MW 35.8 kDa

Uniprot ID **P08174**

NCBI **NP_000565.1**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
CD55
protein