

Product Datasheet

Human B7-2 / CD86 Protein (orb257284)

Description

Human B7-2, His Tag (CD6-H5223) is expressed from



Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target B7-2

Preservatives PBS, pH7.4


Form/Appearance Powder

Storage -20°C


Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.2 kDa. The protein migrates as 40-63 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.



SDS-PAGE analysis of Human B7-2 protein



Line graph illustrates about the sensitivity...

Protein Sequence AAH40261

Purity ● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source Human B7-2, His Tag (orb257284) is expressed from human 293 cells (HEK293). It contains AA Leu 26 - Pro 247 (Accession # AAH40261).

MW 26.2 kDa

Uniprot ID [P42081](#)

NCBI [AAH40261](#)

Expiration Date 6 months from date of receipt.