


Product Datasheet

Human ALCAM Protein (orb257286)

Description	Human ALCAM, His Tag (CD6-H5228) is expressed from human 293 cells (HEK293). It contains AA Trp 28 -...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ALCAM
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 57.8 kDa. The protein migrates as 70-100 kDa under reducing (R) condition, and 66-100 kDa under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001618.2
Purity	● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Source	Human ALCAM, His Tag (orb257286) is expressed from human 293 cells (HEK293). It contains AA Trp 28 - Ala 526 (Accession # Q13740-1).
MW	57.8 kDa
Uniprot ID	Q13740
NCBI	NP_001618.2
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
ALCAM
protein