

Product Datasheet

Human CD27 / TNFRSF7 Protein (orb257294)

Description Human CD27, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Thr 21 - Ile 192 (Acce...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target CD27

Preservatives 50 mM Tris, 100 mM Glycine, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 45.9 kDa. The protein migrates as 53-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_001233

Purity 95%


Source Human CD27, Fc Tag (orb257294) is expressed from human 293 cells (HEK293). It contains AA Thr 21 - Ile 192 (Accession # NP_001233).

MW 45.9 kDa

Uniprot ID **P26842**

NCBI **NP_001233**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
CD27
protein