



Product Datasheet

Human CD47 Protein (orb257295)

Description	Leukocyte surface antigen CD47 is also known as Antigenic surface determinant protein OA3, Integrin-...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD47
Preservatives	50 mM Tris, 100 mM Glycine, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 40.4 kDa. The protein migrates as 55-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_942088
Purity	95%
Source	Human CD47, Fc Tag (orb257295) is expressed from human 293 cells (HEK293). It contains AA Gln 19 - Pro 139 (Accession # NP_942088).
MW	40.4 kDa
NCBI	NP_942088
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human CD47 protein



Line graph illustrates about the sensitivity...