


Product Datasheet

Human CEACAM-1 / CD66a Protein (orb257335)

Description	Human CEACAM-1, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 35 - Gly 428 ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CEACAM-1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 44.1 kDa. The protein migrates as 57-100 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to different glycosylation.
Protein Sequence	AAH14473
Purity	92%
Source	Human CEACAM-1, His Tag (orb257335) is expressed from human 293 cells (HEK293). It contains AA Gln 35 - Gly 428 (Accession # AAH14473).
MW	44.1 kDa
Uniprot ID	P13688
NCBI	AAH14473
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
CEACAM-1
prot...