


Product Datasheet

Human CEACAM-6 / CD66c Protein (orb257336)

Description	Carcinoembryonic antigen-related cell adhesion molecule 6 (non-specific cross reacting antigen) (CEA...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CEACAM-6
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 32.0 kDa. The protein migrates as 45-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_002474.3
Purity	95%
Source	Human CEACAM-6, His Tag (orb257336) is expressed from human 293 cells (HEK293). It contains AA Lys 35 - Gly 320 (Accession # NP_002474.3).
MW	32.0 kDa
NCBI	NP_002474.3
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
CEACAM-6
prot...