


Product Datasheet

Human CEACAM-8 / CD66b Protein (orb257338)

Description	Human CEACAM-8, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 35 - Ser 319 ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CEACAM-8
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 32.2 kDa. The protein migrates as 53-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001807.2
Purity	95%
Source	Human CEACAM-8, His Tag (orb257338) is expressed from human 293 cells (HEK293). It contains AA Gln 35 - Ser 319 (Accession # NP_001807.2).
MW	32.2 kDa
Uniprot ID	P31997
NCBI	NP_001807.2
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
CEACAM-8
prot...