

Product Datasheet

Human C-Reactive Protein (CRP) Protein (orb257363)

Description C-reactive protein (CRP) is a member of the pentraxin family of proteins that are characterized by a...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target C-Reactive Protein

Preservatives Tris-NaCl buffer pH8.0

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 23.9 kDa. The protein migrates as 26 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence NP_000558.2


Purity 95%

Source Human C-Reactive Protein, His Tag (orb257363) is expressed from human 293 cells (HEK293). It contains AA Gln 19 - Pro 224 (Accession # NP_000558.2).

MW 23.9 kDa

NCBI [NP_000558.2](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human C-Reactive pr...