

Product Datasheet

Human M-CSF R / CSF1R / CD115 Protein (orb257372)

Description

Human M-CSF R Protein, Fc Tag (CSR-H5258) is


Reactivity

Human

Endotoxins

0.1 EU per µg

Conjugation

Unconjugated

Target

M-CSF R

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 80.6 kDa. The protein migrates as 100-150 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human M-CSF R Protein, Fc Tag (orb257372) is expressed from human 293 cells (HEK293). It contains AA Ile 20 - Glu 512 (Accession # NP_005202.2).

MW

80.6 kDa

Uniprot ID


P07333

NCBI


NP_005202.2

Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human M-CSF R prote...



Human M-CSF Protein, His Tag (Cat. No. o...