

Product Datasheet

Human CXCR4 / CD184 Protein (orb257389)

Description C-X-C chemokine receptor type 4 is also known as fusin or CD184 (cluster of differentiation 184), CX...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target CXCR4

Preservatives 50 mM Tris, 100 mM Glycine, pH7.5

Form/Appearance Powder

Storage -20°C

Tag N-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 32.3 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAH20968.1


Purity 90%

Source Human CXCR4, Fc Tag (orb257389) is expressed from human 293 cells (HEK293). It contains AA Met 1 - Ser 46 (Accession # AAH20968.1).

MW 32.3 kDa

NCBI [AAH20968.1](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
CXCR4
protein