

Product Datasheet

Human Decorin Protein (orb257393)

Description Human Decorin, His Tag is expressed from human 293 cells (HEK293). It contains AA Gly 17 - Lys 359 (...)

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Decorin

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 38.8 kDa. The protein migrates as 46-55 kDa under reducing (R) condition (SDS-PAGE) due to different glycosylation.

Protein Sequence NP_001911.1

Purity 90%


Source Human Decorin, His Tag (orb257393) is expressed from human 293 cells (HEK293). It contains AA Gly 17 - Lys 359 (Accession # NP_001911.1).

MW 38.8 kDa

Uniprot ID **P07585**

NCBI **NP_001911.1**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
Decorin
prote...