

Product Datasheet

Human DPPIV / CD26 Protein (orb257402)

Description Dipeptidyl peptidase-IV (DPPIV) is also known as adenosine deaminase complexing protein 2, DPPIV or ...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target DPPIV

Preservatives 50 mM HEPES, 150 mM NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag Native

Note For research use only

Application notes This protein carries no "tag". The protein has a calculated MW of 85.5 kDa. The protein migrates as 90-116 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_001926.2


Purity 95%

Source Human DPPIV, Tag Free (orb257402) is expressed from human 293 cells (HEK293). It contains AA Asn 29 - Pro 766 (Accession # NP_001926.2).

MW 85.5 kDa

NCBI [NP_001926.2](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
DPPIV
protein