

Product Datasheet

Human DPPIV / CD26 Protein (orb257403)

Description Dipeptidyl peptidase-IV (DPPIV) is also known as adenosine deaminase complexing protein 2, DPPIV or ...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target DPPIV

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 85.7 kDa. The protein migrates as 95 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_001926.2


Purity 95%

Source Human DPPIV, His Tag (orb257403) is expressed from human 293 cells (HEK293). It contains AA Asp 34 - Pro 766 (Accession # NP_001926.2).

MW 85.7 kDa

NCBI [NP_001926.2](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
DPPIV
protein