

## Product Datasheet

### Human DPPIV / CD26 Protein (orb257404)

**Description** Dipeptidyl peptidase-IV (DPPIV) is also known as adenosine deaminase complexing protein 2, DPPIV or ...

**Reactivity** Human

**Endotoxins** 1.0 EU per µg

**Conjugation** Unconjugated

**Target** DPPIV

**Preservatives** 50 mM Tris, 100 mM Glycine, pH7.5

**Form/Appearance** Powder

**Storage** -20°C

**Tag** N-Fc

**Note** For research use only

**Application notes** This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 112.3 kDa. The protein migrates as 90-116 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Protein Sequence** NP\_001926.2


**Purity** 95%

**Source** Human DPPIV, Fc Tag (orb257404) is expressed from human 293 cells (HEK293). It contains AA Asn 29 - Pro 766 (Accession # NP\_001926.2).

**MW** 112.3 kDa

**NCBI** [NP\\_001926.2](#)

**Expiration Date** 6 months from date of receipt.



SDS-PAGE  
analysis of  
Human  
DPPIV  
protein