


Product Datasheet

Human Ephrin-A3 / EFNA3 Protein (orb257407)

Description	Human Ephrin-A3, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 23 - Ser 213 ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Ephrin-A3
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 48.4 kDa. The protein migrates as 60-65 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAI10407
Purity	90%
Source	Human Ephrin-A3, Fc Tag (orb257407) is expressed from human 293 cells (HEK293). It contains AA Gln 23 - Ser 213 (Accession # AAI10407).
MW	48.4 kDa
Uniprot ID	P52797
NCBI	AAI10407
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Ephrin-A3 pro...