

Product Datasheet

Human Ephrin-A4 / EFNA4 Protein (orb257409)

Description Human Ephrin-A4, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 26 - Gly 171 ...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Ephrin-A4

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 43.0 kDa. The protein migrates as 45-50 kDa and 66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAI07484

Purity 95%


Source Human Ephrin-A4, Fc Tag (orb257409) is expressed from human 293 cells (HEK293). It contains AA Leu 26 - Gly 171 (Accession # AAI07484).

MW 43.0 kDa

Uniprot ID **P52798**

NCBI **AAI07484**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
Ephrin-A4
pro...