

Product Datasheet


Human E-Cadherin / Cadherin-1 Protein (orb257411)

Description

Human E-Cadherin, Fc Tag is expressed from

kDa + -

	Pro 621...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	E-Cadherin
Preservatives	50 mM Tris, 100 mM Glycine, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 92.3 kDa. The protein migrates as 100 kDa and 120 kDa (Cadherin), and 48 kDa (Propeptide) under reducing (R) condition (SDS-PAGE) due to different glycosylation.
Protein Sequence	AAI41839
Purity	95%
Source	Human E-Cadherin, Fc Tag (orb257411) is expressed from human 293 cells (HEK293). It contains AA Gln 23 - Pro 621 (Accession # AAI41839).
MW	92.3 kDa
Uniprot ID	P12830
NCBI	AAI41839
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human E-Cadherin pr...