

Product Datasheet

Human Ephrin-A1 / EFNA1 Protein (orb257414)

Description Human Ephrin-A1, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Asp 19 - Ser 182 ...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Ephrin-A1

Preservatives 50 mM Tris, 100 mM Glycine, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 46.0 kDa. As a result of glycosylation, the protein migrates as 53-57 kDa under reducing (R) condition, and 100-115 kDa under non-reducing (NR) condition (SDS-PAGE).

Protein Sequence NP_004419.2

Purity 95%


Source Human Ephrin-A1, Fc Tag (orb257414) is expressed from human 293 cells (HEK293). It contains AA Asp 19 - Ser 182 (Accession # P20827-1).

MW 46.0 kDa

Uniprot ID **P20827**

NCBI **NP_004419.2**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
Ephrin-A1
pro...