


Product Datasheet

Mouse EGF R Protein (orb257420)

Description	Mouse EGF R, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Ser 647 (Ac...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	EGF R
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 71.2 kDa. The protein migrates as 86-110 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_997538.1
Purity	95%
Source	Mouse EGF R, His Tag (orb257420) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Ser 647 (Accession # Q01279-1).
MW	71.2 kDa
Uniprot ID	Q01279
NCBI	NP_997538.1
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Mouse EGF R
protein