


Product Datasheet

Human EphB4 Protein (orb257423)

Description	Human EphB4, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 16 - Ala 539 (Ac...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	EphB4
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 57.9 kDa. The protein migrates as 65-75 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_004435.3
Purity	97%
Source	Human EphB4, His Tag (orb257423) is expressed from human 293 cells (HEK293). It contains AA Leu 16 - Ala 539 (Accession # NP_004435.3).
MW	57.9 kDa
Uniprot ID	P54760
NCBI	NP_004435.3
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
EphB4
protein