

## Product Datasheet

### Human FABP1 / L-FABP Protein (orb257446)

**Description** Human FABP1, His Tag is expressed from E. coli cells. It contains AA Ser 2 - Ile 127 (Accession # NP...

**Reactivity** Human

**Endotoxins** 1.0 EU per µg

**Conjugation** Unconjugated

**Target** FABP1

**Preservatives** PBS, pH7.4

**Form/Appearance** Powder

**Storage** -20°C

**Tag** N-6×His

**Note** For research use only

**Application notes** This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 15.0 kDa. The protein migrates as 15-16 kDa under reducing (R) condition (SDS-PAGE).

**Protein Sequence** NP\_001434

**Purity** 95%


**Source** Human FABP1, His Tag (orb257446) is expressed from E. coli cells. It contains AA Ser 2 - Ile 127 (Accession # NP\_001434).

**MW** 15.0 kDa

**Uniprot ID** **P07148**

**NCBI** **NP\_001434**

**Expiration Date** 6 months from date of receipt.



SDS-PAGE  
analysis of  
Human  
FABP1  
protein