

Product Datasheet

Human FABP2 / I-FABP Protein (orb257447)

Description Human FABP2, His Tag is expressed from E. coli cells. It contains AA Ala 2 - Asp 132 (Accession # P1...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target FABP2

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag N-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 16.1 kDa. The protein migrates as 15 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence NP_000125.2

Purity 98%


Source Human FABP2, His Tag (orb257447) is expressed from E. coli cells. It contains AA Ala 2 - Asp 132 (Accession # P12104-1).

MW 16.1 kDa

Uniprot ID **P12104**

NCBI **NP_000125.2**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
FABP2
protein