


Product Datasheet

Human FABP3 / H-FABP Protein (orb257448)

Description	Fatty acid-binding protein 3 (FABP3) is also known as Heart-type fatty acid binding protein (H-FABP)...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	FABP3
Preservatives	50 mM Tris, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 15.7 kDa. The protein migrates as 15-17 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_004093
Purity	95%
Source	Human FABP3, His Tag (orb257448) is expressed from E. coli cells. It contains AA Met 1 - Ala 133 (Accession # NP_004093).
MW	15.7 kDa
NCBI	NP_004093
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
FABP3
protein