

Product Datasheet

Human FABP5 / E-FABP Protein (orb257450)

Description Human FABP5, His Tag is expressed from E. coli cells. It contains AA Ala 2 - Glu 135 (Accession # AA...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target FABP5

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag N-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 16.0 kDa. The protein migrates as 16 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence AAH19385

Purity 95%


Source Human FABP5, His Tag (orb257450) is expressed from E. coli cells. It contains AA Ala 2 - Glu 135 (Accession # AAH19385).

MW 16.0 kDa

Uniprot ID **Q01469**

NCBI **AAH19385**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
FABP5
protein