


Product Datasheet

Human Fc Protein (orb257462)

Description	Crystallizable fragments composed of the carboxy-terminal halves of both IMMUNOGLOBULIN HEAVY CHAINS...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IgG1 Fc
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 26.1 kDa. The protein migrates as 32-34 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Glu 99 - Lys 330
Purity	95%
Source	Human IgG Fc, Tag Free (orb257462) is expressed from human 293 cells (HEK293). It contains AA Glu 99 - Lys 330 (Accession # P01857-1).
MW	26.1 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human IgG
Fc protei...