

Product Datasheet

Human Fetuin A / FETUA Protein (orb257467)

Description Alpha-2-HS-glycoprotein (AHSG) is also known as Alpha-2-Z-globulin, Ba-alpha-2-glycoprotein, Fetuin-...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Fetuin A

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 38.2 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_001613.2


Purity 98%

Source Human Fetuin A, His Tag (orb257467) is expressed from human 293 cells (HEK293). It contains AA Ala 19 - Val 367 (Accession # P02765-1).

MW 38.2 kDa

NCBI [NP_001613.2](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
Fetuin A
prot...