


Product Datasheet

Human FGF R4 / CD334 Protein (orb257472)

Description

Human FGF R4, His Tag (FG4-H5228) is expressed


Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

FGF R4

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-6×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 39.3 kDa. The protein migrates as 50-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_002002.3

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source


Human FGF R4, His Tag (orb257472) is expressed from human 293 cells (HEK293). It contains AA Leu 22 - Asp 369 (Accession # P22455-1).

MW

39.3 kDa

Uniprot ID
P22455
NCBI
NP_002002.3
Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human FGF R4 protei...