


Product Datasheet

Human Galectin-3 Protein (orb257492)

Description	Human Galectin-3, His Tag is expressed from E. coli cells. It contains AA Met 1 - Ile 250 (Accession...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Galectin-3
Preservatives	50 mM HEPES, 150 mM NaCl, 5 mM EDTA, 2 mM TCEP, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 27.0 kDa. The protein migrates as 32 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_002297
Purity	92%
Source	Human Galectin-3, His Tag (orb257492) is expressed from E.coli cells. It contains AA Met 1 - Ile 250 (Accession # NP_002297).
MW	27.0 kDa
Uniprot ID	P17931
NCBI	NP_002297
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
Galectin-3
pr...