



Product Datasheet

Human G-CSF R / CD114 Protein (orb257498)

Description	Human G-CSF R, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Glu 25 - Pro 621 (A...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	G-CSF R
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 92.7 kDa. The protein migrates as 100-120 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human G-CSF R, Fc Tag (orb257498) is expressed from human 293 cells (HEK293). It contains AA Glu 25 - Pro 621 (Accession # NP_000751.1).
MW	92.7 kDa
Uniprot ID	Q99062
NCBI	NP_000751.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human G-CSF R protein...



Immobilized Human G-CSF Protein, premium...