



Product Datasheet

SMNDC1 antibody (orb25755)

Description	SMNDC1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP
Immunogen	Recombinant Human Survival of motor neuron-related-splicing factor 30 protein (1-238AA)
Target	SMNDC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075940
Dilution Range	IHC-P:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human he...



Immunoprecipitation analysis of HeLa who...