

Product Datasheet

Human SOST / Sclerostin Protein, His Tag (orb257550)

Description Human SOST, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Tyr 213 (Acc...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target SOST

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag N-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 22.3 kDa. The protein migrates as 30-34 kDa under reducing (R) condition (SDS-PAGE) due to different glycosylation.

Protein Sequence AAI01087.1

Purity 90%


Source Human SOST, His Tag (orb257550) is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Tyr 213 (Accession # AAI01087.1).

MW 22.3 kDa

Uniprot ID [Q9BQB4](#)

NCBI [AAI01087.1](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
SOST
protein