

Product Datasheet

Human HVEM / TNFRSF14 Protein (orb257551)

Description Human HVEM, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 39 - Val 202 (Acce...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target HVEM

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 43.5 kDa. The protein migrates as 55-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_003811.2

Purity 95%


Source Human HVEM, Fc Tag (orb257551) is expressed from human 293 cells (HEK293). It contains AA Leu 39 - Val 202 (Accession # NP_003811.2).

MW 43.5 kDa

Uniprot ID [Q92956](#)

NCBI [NP_003811.2](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
HVEM
protein